



Philips Sonicare
HealthyWhite+
Rechargeable sonic
toothbrush

2 modes

6 customisable settings
1 brush head

HealthyWhite+



HX8911/04

Removes up to 100% more stains in only 1 week*

*than a manual toothbrush

Safely enhance your smile every day without giving up the staining foods and drinks you love! The Philips Sonicare Whitening rechargeable toothbrush is proven to remove up to 100% more stains for whiter teeth in only 1 week*.

Proven to improve oral health

- Removes up to 100% more stains than a manual toothbrush*
- Removes up to 7X more plaque than a manual toothbrush
- Improves gum health in only two weeks

Offers a customised brushing experience

- 6 customisable settings
- Clean and White Mode: proven to remove stains

Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

Provides a superior clean

- DiamondClean brush head offering the best whitening
- Sonicare dynamic cleaning action drives fluid between teeth


PHILIPS
sonicare

Rechargeable sonic toothbrush
2 modes 6 customisable settings, 1 brush head

HX8911/04


Highlights

Whitens teeth




This Philips Sonicare toothbrush removes up to 100% more stains for whiter teeth in just 1 week.

DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond shaped bristles that are medium-stiff to effectively yet gently scrape away plaque. This Philips Sonicare electric toothbrush provides an exceptional cleaning and whiter teeth compared to a manual toothbrush.

Clean and White mode



2-minute Clean mode offers the ultimate in plaque removal (Dental professional recommended brushing time). White mode removes surface stains, brightens and polishes teeth. Removes everyday stains, such as coffee,

tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.

6 customisable settings



Philips Sonicare Stain Removal toothbrush allows you to mix and match modes (Clean mode, White mode) and 3 intensities for maximum comfort and your best clean.

Removes up to 7X more plaque



Removes up to 7 times more plaque from hard-to-reach places than a manual toothbrush.


Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks. Giving a superior cleaning in

between teeth and removes significantly more plaque than a manual toothbrush

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Smartimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Patented sonic technology



The Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Specifications

Modes

- Clean: For exceptional everyday cleaning
- White: Removes surface stains
- 3 intensities: Low, Medium, High

Items included

- Handles: 1 HealthyWhite+
- Brush heads: 1 DiamondClean standard
- Charger: 1

Design and finishing

- Colour: Frost White

Cleaning performance

- Speed: Up to 31,000 brush strokes/min
- Performance: Removes up to 100% more plaque
- Health benefits: Improves gum health in only two weeks
- Whitening benefits: Whitens up to 2 shades in 2

- weeks
- Timer: Smartimer and Quadpacer

Ease of use

- Brush head system: Easy snap-on brush heads
- Brushing time: Up to 3 weeks**
- Battery indicator: Light shows battery status
- Handle: Slim ergonomic design

Technical specifications

- Battery: Rechargeable
- Operating time (full to empty): Up to 3 weeks
- Battery type: Lithium ION

Power

- Voltage: 110-220 V

Service

- Warranty: 2-year limited warranty



Issue date 2015-07-10

Version: 3.1.1

EAN: 08 71010 37145 07

© 2015 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* than a manual toothbrush

**based on two periods of two minutes of brushing per day, in clean mode