



Philips Lumea
IPL hair removal system

Precision Plus

For treatment on body and face
100,000 flashes

Lumea



SC2003/00

Prevent hair regrowth

IPL hair removal: enjoy smooth skin every day

The Philips Lumea IPL hair removal system works wonders to prevent the reappearance of hair on the face as well as the body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

Long lasting smoothness

- Enjoy smooth skin every day
- Effortlessly effective

Hair removal with IPL technology

- Adapted for safe use in the comfort of your own home
- Suitable for a wide variety of hair and skin types

Gentle and effective application for face and body

- Five adjustable light energy settings
- Safe and effective even on sensitive areas
- Precision attachment for safe facial treatment
- Bigger body attachment for fast treatment
- Slide and Flash mode for more convenient application

Maximum freedom, minimum hassle

- Cordless convenience
- No replacement parts, no hidden costs



asimpleswitch.com

PHILIPS

Highlights

Effortlessly effective



Our clinical studies have shown significant hair reduction in as few as four bi-weekly treatments, resulting in smooth, hair-free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth.

IPL technology



Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 women.

Suitable hair and skin types



Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blonde or white/grey hair. Philips Lumea is also NOT suitable for very dark skin.

Five light energy settings



Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.

Application areas




Use with confidence to remove hair on face (upper lip, chin and sideburns) and body parts including legs, underarm, bikini area, stomach and arms.

Cordless convenience



Lumea's unique cordless design is perfect for maximum reach and ease of use.

Extra safety for the face



The precision attachment is equipped with an extra integrated light filter and enables a safe and convenient application on the upper lip, chin and sideburns. It can also be used on other sensitive areas on the body, such as the underarms, the bikini area and hard-to-reach areas.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Technical specifications

- Light fluency (J/cm²): Body attachment: 2-6 J/cm², Precision attachment: 2-6 J/cm
- Light spectrum: Body attachment > 570 nm, Precision attachment > 600 nm
- High performance lamp: Delivers > 100,000 flashes

Safety and adjustable settings

- Safety ring system: Prevent unintentional flashing
- 5 light energy settings: Adjustable to your skin type
- Integrated UV filter: Protects skin from UV light

Application areas

- Body areas: Legs, Underarms, Belly, Bikini
- Face areas: Chin, Sideburns, Upper lip

Attachments

- Body attachment (4 cm²): For use below the neckline
- Precision attachment (2 cm²): For safe application on face

Application time

- Bikini: 4 min
- Face: 1 min
- Lower leg: 8 min
- Slide and Flash: For convenient application
- Step and Flash: For treatment on curvy areas
- Underarms: 1 min

Items included

- Storage: Beauty case
- Instructions for use: Full user manual, Quick Start Guide
- Charger: Multi voltage charger

Power

- Power system: Rechargeable battery
- Voltage: 100-240 V
- Charging time: 100 minutes
- Running time: Minimum 370 flashes

Service

- Warranty: 2 years global warranty



Issue date 2014-01-16

Version: 7.1.1

EAN: 08 71010 35794 65

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com