



Philips Sonicare  
DiamondClean  
Standard sonic toothbrush  
heads

4-pack  
Standard

**DiamondClean**

**HX6064/33**



# Superior cleaning. Whiter teeth.

Toothbrush head for whiter teeth

Replacement brush head for Sonicare electric toothbrush. Best brush head for an ultimate clean

#### **Sonicare's best plaque removal**

- Removes up to 100% more plaque than a manual toothbrush\*

#### **Provides a superior clean**

- Diamond-Shaped bristles sweep away plaque
- 44% more bristles for a more thorough clean
- Medium stiff bristles provide a firm, but gentle experience.

#### **Designed to optimise performance**

- Reminder bristles ensure you are getting an effective clean

#### **A brush head that fits multiple handles**

- Snap-on brush head

**PHILIPS**  
**sonicare**

# Specifications

## Ease of use

- Suitable for these models: FlexCare+, FlexCare, HealthyWhite, HydroClean, EasyClean, DiamondClean
- Brush head system: Easy snap-on brush heads for

optimal hygiene

## Technical specifications

- Operating time: For optimal results, a new brush head is recommended every three months.

# Highlights

## 44% more bristles

44% more bristles for a more thorough clean.

## Diamond-Shaped bristles

Diamond-Shaped bristles create more scraping surfaces to sweep away plaque.

## Medium stiff bristle

Medium stiff bristles provide a firm, but gentle experience.


## Reminder bristles

Reminder bristles let you know when to replace the brush head. After three months of normal use bristles exhibit fatigue, and brush heads are less effective. Replace your brush head every 3 months.

## Removes up to 100% more plaque

\*Removes up to 100% more plaque from hard-to-reach places than a manual toothbrush

## Snap-on brush head



Easy brush head placement and handle cleaning



Issue date 2013-12-09

Version: 1.0.1

EAN: 08 71010 36549 19

© 2013 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)