



Philips Sonicare Sensitive  
Standard sonic toothbrush  
heads

2-pack



HX6052/26

## Gentle cleaning

### Gentle Sonicare toothbrush head

Replacement brush head for the Sonicare electric toothbrush. Soft on teeth and gums. Effective on plaque.

#### **Softest Sonicare brush head ever**

- Innovative new brush head profile delivers a gentle experience
- Ultra-soft bristles for a more effective, gentle clean

#### **Designed to optimise performance**

- Reminder bristles ensure you are getting an effective clean

#### **A brush head that fits multiple handles**

- Snap-on brush head

#### **Gentle yet effective cleaning**

- Gently removes plaque for an absolute clean and smooth feel

**PHILIPS**  
**sonicare**

Standard sonic toothbrush heads  
2-pack

# Specifications

## Ease of use

- Suitable for these models: FlexCare+, FlexCare, HealthyWhite, DiamondClean, EasyClean, HydroClean
- Brush head system: Easy snap-on brush heads for

optimal hygiene

## Technical specifications

- Operating time: For optimal results, a new brush head is recommended every three months.

HX6052/26

# Highlights

## Snap-on brush head



Easy brush head placement and handle cleaning

## Reminder bristles

Reminder bristles let you know when to replace the brush head. After three months of normal use bristles exhibit fatigue, and brush heads are less effective. Replace your brush head every 3 months.

## Gently removes plaque

Gently removes plaque for an absolute clean and smooth feel.

## Innovative new brush head

Innovative new brush head profile delivers a gentle experience while brushing sensitive teeth and gums.

## Ultra-soft bristles

Ultra soft bristles at the core of the bristle field are 25% softer than our standard ProResults brush head. For a more effective, gentle clean.



Issue date 2013-12-09

Version: 1.0.1

EAN: 08 71010 36418 34

© 2013 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)