



Philips Sonicare
DiamondClean
Rechargeable sonic
toothbrush

5 modes

1 brush head
Glass charger
USB travel charger

DiamondClean



HX9331/04

Whiter, healthier teeth for life

Best Sonicare whitening electric toothbrush

Sonicare's best whitening in our most elegant Philips Sonicare electric toothbrush. Switch to Sonicare.

Proven to improve oral health

- Removes up to 7X more plaque than a manual toothbrush
- Improves gum health in only two weeks
- Whitens teeth more than 2x better than a manual toothbrush
- Sensitive Mode: Gentle teeth and gum cleaning
- Philips Sonicare toothbrush helps whiten teeth

Provides a superior clean

- DiamondClean toothbrush head for Sonicare's best whitening
- Sonicare dynamic cleaning action drives fluid between teeth

Designed around you

- Dual Charging System: Charging glass and travel case.

Guides you to follow dental recommendations

- Quadpacer interval timer encourages thorough brushing
- Two-minute timer helps ensure recommended brushing time

PHILIPS

sonicare

Highlights

Removes up to 7X more plaque




Removes up to 100% more plaque from hard to reach places than a manual toothbrush.

Improves gum health



This Philips Sonicare electric toothbrush provides optimal cleaning between teeth and along the gumline for improved gum health in just two weeks. Giving a superior cleaning in between teeth and removes significantly more plaque than a manual toothbrush

Whitens teeth twice as well



This Philips Sonicare toothbrush whitens teeth more than 2x better than a manual toothbrush in just 1 week. Removes up to 100% more stains for whiter teeth in just 1 week.

DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond-shaped, medium stiffness bristles to effectively yet gently scrape away plaque. This Philips Sonicare electric toothbrush provides exceptional cleaning and whiter teeth compared to a manual toothbrush.

Dual Charging System




1) Charger glass - Simply place your Philips Sonicare electric toothbrush in the glass to charge. Can also be used to rinse your mouth after brushing. 2) Charging travel case - Connects to your laptop via USB or can be charged by wall socket. Allows you to take your DiamondClean anywhere and everywhere.

Gentle teeth and gum cleaning



Philips Sonicare toothbrush with Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Smartimer



2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Helps whiten teeth



Philips Sonicare's best teeth whitening: this toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Specifications

Modes

- Clean: For exceptional everyday cleaning
- Gum Care: Gently massages gums
- Polish: Brighten and polish your teeth
- Sensitive: Gentle teeth and gum cleaning
- White: Removes surface stains

Items included

- Handles: 1 DiamondClean
- Brush heads: 1 DiamondClean standard
- Glass charger: 1
- Travel case: USB travel charger

Design and finishing

- Colour: Ceramic white

Cleaning performance

- Speed: Up to 31,000 brush strokes/min
- Performance: Removes up to 100% more plaque
- Health benefits: Improves gum health in only two weeks

- Whitening benefits: Whitens teeth twice as well
- Timer: Smartimer and Quadpacer

Ease of use

- Brush head system: Easy snap-on brush heads
- Brushing time: Up to 3 weeks**
- Battery indicator: Illuminated icon indicates battery life
- Handle: Slim ergonomic design
- Display: Illuminated display

Technical specifications

- Battery: Rechargeable
- Battery type: Lithium ION
- Operating time (full to empty): Up to 3 weeks

Power

- Voltage: 110-220 V

Service

- Warranty: 2-year limited warranty



Issue date 2015-02-27

Version: 2.0.1

EAN: 08 71010 37102 19

© 2015 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

**based on two periods of two minutes of brushing per day, in clean mode