



Philips Sonicare FlexCare  
Rechargeable sonic  
toothbrush

**FlexCare**



**HX6982**

## Advanced Cleaning

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare that adapts to your oral care needs.

### **Guides you to follow dental recommendations**

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

### **Provides a superior clean**

- Dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

### **Proven to improve oral health**

- Naturally whiter teeth with patented sonic technology
- Clinically proven safe and gentle
- Improves gum health in just two weeks

### **Offers you a customised cleaning experience**

- Easy-start programme to acclimate to the Sonicare experience

**PHILIPS**  
**sonicare**

sense and simplicity

# Specifications

## Cleaning performance

- Speed: Up to 31,000 brush strokes per minute

## Items included

- Additional brush head colour code rings: 3 pcs
- Deluxe Charger with cord wrap: 1
- Deluxe soft travel case: 1
- Educational DVD: 1
- FlexCare handle: 1

## Technical specifications

- Battery: Rechargeable
- Operating time (full to empty): Fully charged should provide 20 2-min. brushings

## Design and finishing

- Colour: Medical Blue
- Chrome front panel

## Features

- Deluxe recharge gauge: Light signal helps track battery charge
- Easy-Start: Gently increases power over 14 uses
- Quadpacer: Encourages attentive brushing of each quadrant of the mouth
- Rechargeable battery
- Smartimer: 2-minute timer helps ensure recommended brushing

# Highlights

## Patented sonic technology

Sonicare's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

## Angled brush head neck

Uniquely angled brush head neck makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Naturally whiter teeth

Dynamic Cleaning Action and direct broader contact with each tooth have been proven to remove everyday stains and help keep teeth naturally whiter.

## Safe and gentle

Your Sonicare is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

## Improves gum health

Optimal cleaning between teeth and along the gumline provides improved gum health in just two weeks. Giving a superior clean in between teeth and removes significantly more plaque than a manual toothbrush

## Smartimer

Two-minute timer automatically tracks brushing time and turns the toothbrush off after two full minutes of brushing is complete.

## Quadpacer

30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

## Easy-start programme

Gently increases power over the first 14 uses to ease into the Sonicare experience

