

Philips Lumea
IPL hair removal system

For use on body only

15 minutes to treat lower legs
Lifetime >100,000 light pulses
Cordless design

Lumea



SC2004/11

The ideal way to prevent hair regrowth

for smooth skin at home

Philips Lumea IPL is the most effective cordless IPL solution and works wonders to prevent the reappearance of hair on the body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

Long lasting smoothness

- Enjoy smooth skin every day
- Effortlessly effective

Hair removal with IPL technology

- Adapted for safe use in the comfort of your own home
- Suitable for a wide variety of hair and skin types

Gentle and effective, even on sensitive skin

- Safe and effective even on sensitive areas
- Bigger body attachment for fast treatment
- Five adjustable light energy settings
- Slide and Flash mode for more convenient application

Maximum freedom, minimum hassle

- Cordless convenience
- No replacement parts, no hidden costs



asimpleswitch.com

PHILIPS

IPL hair removal system

For use on body only 15 minutes to treat lower legs, Lifetime >100,000 light pulses, Cordless design

Highlights

Long lasting smoothness



Philips Lumea applies gentle pulses of light to the hair root. This results in the hair shedding naturally and hair regrowth is inhibited. Repeating this treatment regularly will keep your skin smooth every day.

Effortlessly effective



Our clinical studies have shown significant hair reduction in as few as four bi-weekly treatments, resulting in smooth, hair-free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth.

IPL technology



Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop

our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 women.

Suitable hair and skin types

Suitable hair and skin tones

Skin Tone	White	Beige	Light brown	Medium brown	Dark brown	Brownish black and darker
Hair Color						
Black	⊙	⊙	⊙	⊙	⊙	⊙
Dark brown	⊙	⊙	⊙	⊙	⊙	⊙
Brown	⊙	⊙	⊙	⊙	⊙	⊙
Dark blond	⊙	⊙	⊙	⊙	⊙	⊙
Light blond	⊙	⊙	⊙	⊙	⊙	⊙
Red	⊙	⊙	⊙	⊙	⊙	⊙
White/grey	⊙	⊙	⊙	⊙	⊙	⊙

⊙ This product is suitable for you to use ⊙ This product is neither suitable nor safe for your skin and hair type

Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blond or white/grey hair. Philips Lumea is also NOT suitable for very dark skin.

Five light energy settings



Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.


Application areas



Use with confidence to remove hair on body areas below the neckline: legs, underarms, bikini area, arms and belly.

SC2004/11

Cordless convenience



Lumea's unique cordless design is perfect for maximum reach and ease of use.

No replacement parts needed



Philips Lumea is a complete solution that comes ready to work straight out of the box. More importantly, it does not require any replacement lamps or gels.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

IPL hair removal system

For use on body only 15 minutes to treat lower legs, Lifetime >100,000 light pulses, Cordless design

SC2004/11

Specifications

Technical specifications

- Light fluency (J/cm²): 2.5-5.6
- Light spectrum: >570 nm
- High performance lamp: Delivers > 100,000 flashes

Safety and adjustable settings

- 5 light energy settings: Adjustable to your skin type
- Integrated UV filter: Protects skin from UV light
- Integrated safety system: Prevent unintentional flashing

Application areas

- Body areas: Legs, Belly, Bikini, Underarms

Application mode

- Slide and Flash: For convenient application
- Step and Flash: For treatment on curvy areas
- Corded / cordless use: Cordless

Attachments

- Body attachment (4 cm²): For use below the

neckline

Application time

- Lower leg: 8 min
- Underarm: 1 min
- Bikini line: 1 min

Items included

- Storage: Pouch
- Instructions for use: Quick Start Guide, User manual
- Charger: Multi voltage charger

Power

- Voltage: 100-240 V
- Charging time: 100 min
- Running time: Minimum 370 flashes

Service

- Warranty: 2 years global warranty



Issue date 2015-06-23

Version: 7.2.1

EAN: 08 71010 36578 35

© 2015 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com