

A woman with blonde hair tied back, wearing a white turtleneck and a beige knit cardigan, is sitting and holding a baby. The baby is wearing an orange shirt and looking towards the camera. A purple pillow is visible in the background.

Philips Avent
Comfort Single electric
breast pump

Natural
Includes 4 oz bottle



SCF332/01

More comfort, more milk naturally

Breast pump with massage cushion

When you are comfortable and relaxed, your milk flows more easily. That's why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forwards and let our soft massage cushion gently stimulate your milk flow.

Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

Choose your most comfortable setting

- Features a gentle stimulation mode and 3 pumping settings

Gently stimulates natural let down and milk flow

- Soft massage cushion with warm feel

Easy to combine breast and bottle feeding

- Includes our Natural bottle and teat for natural latch on

Other benefits

- Compact lightweight design
- Intuitive assembly. Easy visual matching of parts
- Compatible with other Philips Avent feeding products
- Easy cleaning due to the small number of separate parts
- Easy operation at the touch of a button

Highlights

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping; no need for you to lean forwards to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping naturally helps your milk to flow more easily.

Simple settings to choose from



When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortably for you.

Soft massage cushion



Our massage cushion has a new soft velvety texture that gives a warm feel to the skin for

comfortable, gentle stimulation of your milk flow. During pumping, it gently mimics your baby's suckling for natural let down and clinically proven effectiveness. The soft massage cushion is designed to comfortably fit the majority of mums. Should you need it, a larger size cushion is available separately (availability may vary per country).

Includes our Natural teat



The wide breast-shaped teat promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Compact design



The breast pump has a compact design, which makes it easy to hold and position on your breast. The small, lightweight base unit can easily be placed within comfortable reach for full control when pumping. For extra convenience in transport and storage, the tube simply wraps around the base unit. Easy to use on the go with batteries.

Easy cleaning

Cleaning is easy thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts

are dishwasher-proof, except the electrical parts.

Easy operation

A single electric breast pump is ideal for mums who pump milk frequently and value being able to pump effortlessly, at the touch of a button.

Intuitive assembly



Easy visual matching of parts for intuitive assembly.



[asimpleswitch.com](https://www.asimpleswitch.com)

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Country of origin

- England

Material

- Breast pump: Polypropylene, BPA free*
- Bottle: Polypropylene, BPA free*
- Teat: Silicone, BPA free*

What is included

- Breast pump body: 1 pcs
- Base unit incl. tubing: 1 pcs
- Standard-size cushion: 1 pcs
- Natural bottle 4 oz: 1 pcs
- Extra soft Newborn Flow Teat: 1 pcs
- Travel cover: 1 pcs
- Sealing disc for milk storage: 1 pcs
- Breast pad sample packs: 2 (2 Day pads and 2 Night pads) pcs
- Spare diaphragm: 1 pcs

Breast pump

- Material: BPA free*

Design

- Breast pump design: Compact design
- Bottle design: Ergonomic shape, Wide neck

Ease of use

- Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Development stages

- Stages: Pregnancy, 0–6 months

Functions

- No leaning forwards: Sit in a comfortable position
- Settings: 3 Expression settings, 1 Stimulation mode
- Soft massage cushion: Gentle stimulation



Issue date 2015-06-08

Version: 2.1.1

EAN: 08 71010 35790 83

© 2015 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Following EU directive (2011/8/EU)