

# VH5000E

## Cool Mist – Mini Ultrasonic Humidifier

*For better breathing and a good night 's sleep*



- Help reduce the survival of flu viruses and ease dry nose, throat, eyes and skin discomforts with moisture level between 40%-60% . <sup>1</sup>
- Suitable for babies and children rooms : ultra quiet and no risk of burns
- For more comfort release of Vicks soothing menthol vapors when used with VapoPads® (Optional).
- To help prevent white dust, the water passes through Demineralization Cartridge (Optional).
- For Germ killing function, the water passes through an Anti-Bacterial Cartridge (Optional).
- Change cartridges every 30-60 days, based on local water supply
- Mist output regulation switch to adjust the moisture level



### -----Consumables-----



VapoPads (5 pcs)  
VH5A or VH5B



Demineralisation Cartridge  
CH5000-DEM  
EAN: 4022167501653



Antibacterial Cartridge  
CH5000-ABA  
EAN: 4022167502254

**2** Year  
Guarantee

Studies have shown that keeping indoor air at an RH level of 40-60% reduces survival of flu viruses on surfaces and in the air compared to lower RH levels. <sup>1</sup>

### Technical Information VH5000E1

Article	EAN	4022167500175
Master box	ITF	14022167500714
Input power	Watts	21
Voltage/Frequency	V/~Hz	220-240/~50
Sound level	db(A)	Max 35
Output Performance	l/day	3.7
Tank capacity	l	1.8
Hrs between fillings	(high setting)	8
Ideal for rooms	approx. m²	15
Length of supply cord incl. plug	approx. m	1.83
Color		White
Dimensions of product	L x W x H	23.6 x 16.7 x 25.7 cm
Dimensions of single box	L x W x H	26.3 x 19.3 x 28.5 cm
Weight of single box	approx. kg	2
Master box quantity	pcs.	2
Dimensions master box	L x W x H	40.2 x 27.5 x 31.1
Weight master box	approx. kg	4.7



Specifications are subject to change without notice  
VH5000E

VICKS is a registered trademark of the Procter & Gamble Company  
Manufactured by Kaz Europe SA under license from the Procter & Gamble Company, Cincinnati, Ohio, USA  
© 2009 Kaz Europe SA

**Kaz Europe S.A**  
Healthcare & Home Comfort Products  
Place Chauderon 18  
1003 Lausanne - Switzerland

Phone +41 21 644 01 10  
Fax +41 21 644 01 11  
e-Mail [info-europe@kaz.com](mailto:info-europe@kaz.com)  
Web [www.kaz.com](http://www.kaz.com) [www.kaz-europe.com](http://www.kaz-europe.com)