

IRT3030

Braun ThermoScan® 3

Compact infrared ear thermometer for the whole family


New

- Reliable & easy to use entry ear thermometer
- From no. 1 brand among doctors*
- With **new** design and **new** Audio fever indicator



Fast & accurate

- Precise measurement in just 1 second
- Measures the infrared heat generated by the eardrum and the surrounding tissues, with clinically proven accuracy
- Improved clean ability of the sensor tip if accidentally used without a probe cover (**new**)



New Audio fever indicator

- Easy to interpret the temperature value: number of beeps indicates normal (1 beep), elevated (2 beeps) and high temperature (4 beeps)



Features for easy & pleasant usability

- Compact design is easy to hold
- Easy to read LCD display
- Probe suitable for use with newborns
- Memory recalls the last temperature reading
- Auto shut-off after 60 seconds for longer battery life.
- Readings in Celsius (34 - 42.2° C) or Fahrenheit (93.2 - 108° F)




Disposable lens filters for hygienic use

- With BPA & latex free disposable probe covers i.e. no need to clean or disinfect the probe
- Using a clean disposable cover with each measurement:
 - 1) ensures hygienic use
 - 2) protects from cross-contamination
 - 3) ensures earwax residues do not impact accuracy of the reading.



Accessories

- Probe cover protection cap
- 21x lens filter
- 1 x cell battery. Battery door with a safety lock to prevent choking hazard (**new**)



Impactful in-store presentation

- Attractive packaging with premium look and feel, with visible product. Hanging hole for in store merchandising

* Based on survey of paediatricians and general practitioners in UK, France, Germany and the Netherlands. IPSOS, 2012



Helen of Troy

Kaz Europe Sàrl
Health Care, Personal Care & Home Comfort Products
Place Chauderon 18
1003 Lausanne - Switzerland
Phone +41 21 644 01 10
Fax +41 21 644 01 11
e-Mail info-europe@kaz.com
www.hot-europe.com

BRAUN

IRT3030

Braun ThermoScan® 3

Compact infrared ear thermometer for the whole family

New



	Quantity	Weight (kg)	Dimensions H x L x W (mm)	EAN/TUC Code	Oracle code
Unit		0.05	128 x 31 x 28		IRT3030WE
Unit box	1	0.214	196 x 42 x 98	4022167330307	
Master case	6	1.284	135 x 386 x 216	14022167330304	
Unit	1	0.05	128 x 31 x 28		IRT3030EE
Unit box	1	0.214	196 x 42 x 98	4022176230300	
Master case	6	1.284	135 x 386 x 216	14022167230307	
Nb. layers	13				
Nb. cases per layers	10				
Pallet	780				


NEW

Refreshed design & package



NEW

Reliable, accurate & easy to clean sensor



NEW

More focus on consumable benefits



NEW

Safety lock to protect against the choking hazard of the battery for children




Fast measurements (1 sec)



NEW

Audio fever indicator



Stores the last reading



Suitable for newborns



Helen of Troy

2 year guarantee, CE mark

Kaz Europe Sàrl
 Health Care, Personal Care & Home Comfort Products
 Place Chauderon 18
 1003 Lausanne - Switzerland
 Phone +41 21 644 01 10
 Fax +41 21 644 01 11
 e-Mail info-europe@kaz.com
 www.hot-europe.com

BRAUN