



Philips AVENT
Microwave Steam Steriliser

Microwave Steriliser

SCF281/02



Easy sterilising at home and away

Sterilises 4 Philips AVENT bottles in 2 minutes*

The Philips AVENT Microwave Steam Steriliser SCF281/02 has a lightweight and compact design, which makes it ideal for use in and out of the home. Contents remain sterile for up to 24 hours if the lid is unopened.

Ultra fast, easy to use

- Sterilises your bottles in as little as 2 minutes

Compact and lightweight

- Convenient for travel. Fits most microwaves.

Safe, easy handling

- Side grips close the lid securely

Effective sterilisation

- Natural steam sterilisation kills 99.9% of harmful germs

BPA-Free

- This steriliser is made from BPA-Free materials



asimpleswitch.com

PHILIPS
AVENT

Highlights

Convenient for travel



The Philips AVENT Microwave Steriliser has been designed to fit most microwaves on the market. Its small size makes it convenient for travel, ensuring that you always have a sterile bottle when you are going on a short overnight trip or on a longer holiday abroad. Also ideal as an extra steriliser for at the grandparents' house. Dimensions: 166 (H), 280 (W), 280 (L) mm.

Ready in 2 minutes



Just add water, load and place in the microwave for as little as 2 minutes. The exact length of the cycle will depend on the wattage of your microwave: 2 min at 1100-1850 Watt, 4 min at 850-1000 Watt, 6 min at 500-850 Watt.

Clips close the lid securely



The microwave steriliser has clips on the side for extra safety. The clips close the lid securely so hot

water cannot easily escape when you take the steriliser out of the microwave. The side grips have also been designed to stay cooler to help you handle the steriliser safely.

Kills 99.9% of harmful germs



Sterilising protects your baby from particularly harmful bacteria until his immune system is strong enough. Philips AVENT Sterilisers use the hospital method of steam sterilisation, which is quick, easy and effective, with no chemicals involved. Contents stay sterile for 24 hours if the lid is unopened.

BPA-Free



This steriliser is made from BPA-Free materials

Compatible with Philips AVENT

Compatible with all Philips AVENT products

Specifications

Technical specifications

- Sterilisation Time: 2 min at 1100-1850 W, 4 min at 850-1000 W, 6 min at 500-800 W
- Water capacity: 200 ml

Weight and dimensions

- Weight: 740 g
- Dimensions: 166 (H), 280 (W), 280 (L) mm

Country of origin

- England

What is included

- Microwave steam steriliser: 1 pcs
- Tongs: 1 pcs

Development stages

- Stages: 0-6 months



Issue date 2014-03-27

Version: 2.0.1

EAN: 08 71010 35675 30

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* 2 min at 1100-1850 W, 4 min at 850-1000 W, 6 min at 500-800 W.